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<Result Summary>

Your ImmunAge-based age groups for each measured gene

If your calculated age group is older than your chronological age, improve your lifestyle to rejuvenate your immunity by following the **key phrases** below.

Please refer to ImmunAge Handbook for details.

Gene Name	Copies of this gene decrease as we age (Handbook P.9)
CD248	Resistance to infection : Vitamins and moderate exercise. This gene is highly expressed by naive T cells. The proportion of Naive T cells in the entire immune cell population declines as we age, and this is one of the main reasons why elderly people are more susceptible to infectious diseases.
Your Result	30s

Gene Name	Copies of this gene increase as we age (Handbook P.13)
ISM1	Skin and mucosal immunity: Vitamins and moderate exercise. This gene is related to the activities of NK (Natural Killer) cells and NKT (Natural Killer T) cells. Aside from peripheral blood, ISM1 is also highly expressed in immune cells residing in lung, mucosa membrane such as intestine, and skin ¹⁵ .
Your Result	40s

Gene Name	Copies of this gene increase as we age (Handbook P.16)
NEFL	Brain and cognition : Omega 3 fatty acids and polyphenol. Product of this gene is a neurofilament (a fibrous element found in nerve cells). Activity of NEFL gene increases in response to brain injury and nerve damage ²² . The amount of NEFL/neurofilament in cerebrospinal fluid inversely correlates with the score of mental state examination ²⁷ .
Your Result	50s

Gene Name	Copies of this gene increase as we age (Handbook P.19)
CRIP1	Zinc intake: Zinc-rich foods This gene is highly expressed in peritoneal macrophages and in peripheral blood immune cells. Over-activation of CRIP1 may cause immune system dysfunction. Protein product of CRIP1 binds to zinc, and is highly expressed in intestinal mucosa ⁴³ . Zinc deficiency negatively affects the function of immune system ^{41,46} .
Your Result	60s

Gene Name	Copies of this gene decrease as we age (Handbook P.21)
SLC4A10	Salt balance: Water intake and dietary salt This gene encodes for bicarbonate transporter, which exchanges ions across cell membrane and contributes to maintain pH balance between inner and outer of a cell. Also, functional capability of this class of transporter influences the relationship between dietary salt intake and blood pressure ^{22,29} .
Your Result	30s

※Numbers in the text correspond to the reference literatures listed in ImmunAge Handbook.

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