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**1. Transient appearance of circulating tumor DNA associated with de novo treatment.**

Sci Rep. 2016 Dec 9;6:38639. doi: 10.1038/srep38639. PubMed PMID: 27934896; PubMed Central PMCID: PMC5146655.

Kato K, Uchida J, Kukita Y, Kumagai T, Nishino K, Inoue T, Kimura M, Imamura F.

**2. Early responses of EGFR circulating tumor DNA to EGFR tyrosine kinase inhibitors in lung cancer treatment.**

Oncotarget. 2016 Nov 1;7(44):71782-71789. doi: 10.18632/oncotarget.12373.

Imamura F, Uchida J, Kukita Y, Kumagai T, Nishino K, Inoue T, Kimura M, Kato K.

**3. Numerical indices based on circulating tumor DNA for the evaluation of therapeutic response and disease progression in lung cancer patients.** Sci Rep. 2016 Jul 6;6:29093. doi: 10.1038/srep29093

Kato K, Uchida J, Kukita Y, Kumagai T, Nishino K, Inoue T, Kimura M, Oba S, Imamura F

**4. Monitoring of treatment responses and clonal evolution of tumor cells by circulating tumor DNA of heterogeneous mutant EGFR genes in lung cancer.**

Lung Cancer. 2016 Apr;94:68-73. doi: 10.1016/j.lungcan.2016.01.023. Epub 2016 Feb 2

Imamura F, Uchida J, Kukita Y, Kumagai T, Nishino K, Inoue T, Kimura M, Oba S, Kato K.

**5. Dynamics of circulating tumor DNA represented by the activating and resistant mutations in epidermal growth factor receptor tyrosine kinase inhibitor treatment**

Cancer Sci. 2016 Mar;107(3):353-8. doi: 10.1111/cas.12860. Epub 2016 Feb 18

Uchida J, Imamura F, Kukita Y, Oba S, Kumagai T, Nishino K, Inoue T, Kimura M, Kato K.

**6. Diagnostic Accuracy of Noninvasive Genotyping of EGFR in Lung Cancer Patients by Deep Sequencing of Plasma Cell-Free DNA.**

Clin Chem. 2015 Sep;61(9):1191-6. doi: 10.1373/clinchem.2015.241414. Epub 2015

Uchida J, Kato K, Kukita Y, Kumagai T, Nishino K, Daga H, Nagatomo I, Inoue T, Kimura M, Oba S, Ito Y, Takeda K, Imamura F.

**7. Quantitative identification of mutant alleles derived from lung cancer in plasma cell-free DNA via anomaly detection using deep sequencing data.**

PLoS One. 2013 Nov 21;8(11):e81468. doi: 10.1371/journal.pone.0081468.

Kukita Y, Uchida J, Oba S, Nishino K, Kumagai T, Taniguchi K, Okuyama T, Imamura F, Kato K

The 57th Annual Meeting of the Japan Lung Cancer Society

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### <Result Summary>

#### Your ImmunAge-based age groups for each measured gene

If your calculated age group is older than your chronological age, improve your lifestyle to rejuvenate your immunity by following the **key phrases** below.

Please refer to ImmunAge Handbook for details.

Gene Name	Copies of this gene decrease as we age (Handbook P.9)
<b>CD248</b>	<b>Resistance to infection : Vitamins and moderate exercise.</b> This gene is highly expressed by naive T cells. The proportion of Naive T cells in the entire immune cell population declines as we age, and this is one of the main reasons why elderly people are more susceptible to infectious diseases.
Your Result	<b>30s</b>

Gene Name	Copies of this gene increase as we age (Handbook P.13)
<b>ISM1</b>	<b>Skin and mucosal immunity: Vitamins and moderate exercise.</b> This gene is related to the activities of NK (Natural Killer) cells and NKT (Natural Killer T) cells. Aside from peripheral blood, ISM1 is also highly expressed in immune cells residing in lung, mucosa membrane such as intestine, and skin <sup>15</sup> .
Your Result	<b>40s</b>

Gene Name	Copies of this gene increase as we age (Handbook P.16)
<b>NEFL</b>	<b>Brain and cognition : Omega 3 fatty acids and polyphenol.</b> Product of this gene is a neurofilament (a fibrous element found in nerve cells). Activity of NEFL gene increases in response to brain injury and nerve damage <sup>22</sup> . The amount of NEFL/neurofilament in cerebrospinal fluid inversely correlates with the score of mental state examination <sup>27</sup> .
Your Result	<b>50s</b>

Gene Name	Copies of this gene increase as we age (Handbook P.19)
<b>CRIP1</b>	<b>Zinc intake: Zinc-rich foods</b> This gene is highly expressed in peritoneal macrophages and in peripheral blood immune cells. Over-activation of CRIP1 may cause immune system dysfunction. Protein product of CRIP1 binds to zinc, and is highly expressed in intestinal mucosa <sup>43</sup> . Zinc deficiency negatively affects the function of immune system <sup>41,48</sup> .
Your Result	<b>60s</b>

Gene Name	Copies of this gene decrease as we age (Handbook P.21)
<b>SLC4A10</b>	<b>Salt balance: Water intake and dietary salt</b> This gene encodes for bicarbonate transporter, which exchanges ions across cell membrane and contributes to maintain pH balance between inner and outer of a cell. Also, functional capability of this class of transporter influences the relationship between dietary salt intake and blood pressure <sup>52,55</sup> .
Your Result	<b>30s</b>

※Numbers in the text correspond to the reference literatures listed in ImmunAge Handbook.

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