AGENDIA
PRECISION ANSWERS

MammaPrint + BluePrint 検査報告書ガイド

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Summary of Results

PATIENT NAME: LAST NAME, FIRST NAME
DOB: 10-JAN-1961

GENDER: Female
SPECIMEN ID: MRN 123456
PATIENT MRN: 825826362
CUSTOMER REF: 123456789

ORDERED BY: Dr John Doe
ACCOUNT: John Doe Hospital
1234 Main St.
Irvine CA 92618 USA

REQUISITION #: 000
SPECIMEN TYPE: FFPE, Core
SPECIMEN SOURCE: Left Breast
COLLECTED DATE: 18-Feb-2014
RECEIVED DATE: 19-Feb-2014
REPORTED DATE: 21-Feb-2014

Summary of Results: LOW RISK LUMINAL-TYPE (A)

MammaPrint Risk-of-Recurrence Result: **LOW RISK**

MammaPrint Index: +0.455
Low Risk (p < 0.0001)

Average 15-Year Risk-of-Recurrence if Untreated: 10%

Blueprint Molecular Subtype Result: **LUMINAL-TYPE**

Predicted Risk of Recurrence Without Adjuvant Systemic Treatment

The integration of clinical risk assessment with MammaPrint results can help refine an individual's prognosis to help better guide the most appropriate treatment strategy. The percentage of patients without distant recurrence at 5 years (DMFS) shown in the diagrams below were observed in the MINDACT trial (clinical risk can be determined by utilizing the clinical risk algorithm on page 2).

Clinical Pathological Low Risk/MammaPrint Low Risk³

% of patients without metastasis at 5 years (Distant Metastasis Free Survival)

5 Year Endocrine Therapy Alone (if applicable): 97.6%

Clinical Pathological High Risk/MammaPrint Low Risk³

% of patients without metastasis at 5 years (Distant Metastasis Free Survival)

Endocrine Therapy Alone: 94.4%

Endocrine Therapy + Chemotherapy: 95.9%

A MammaPrint Low Risk result does not guarantee that cancer will not recur. A High Risk result does not guarantee that cancer will recur. Individual results may vary. Many other results should be taken in the context of other relevant clinicopathological factors and standard practice of medicine.
*A result is the likelihood of achieving a clinical result by chance, assuming there is no difference between the two arms being investigated.
The test is not intended to be used for the purpose of clinical trial stratification or for the purpose of drug development.

Agenda NV, Science Park 406, 3508 XH Amstelveen, The Netherlands T +31 (0)20 452 1010 customer.service@agenda.com Page 1 of 3

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1 Summary of Results

PATIENT NAME: LAST NAME, FIRST NAME
DOB: 10-JAN-1961

GENDER: Female
SPECIMEN ID: MSN 123456
PATIENTURN: 94292922
CUSTOMER REF: 123456789

ORDERED BY: Dr. John Doe
ACCOUNT: John Doe Hospital
1234 Main St.
Irvine CA 92618 USA

REQUESTION #: 000
SPECIMEN TYPE: FFPE, Core
SPECIMEN SOURCE: Left Breast
COLLECTED DATE: 18-Feb-2014
RECEIVED DATE: 18-Feb-2014
REPORTED DATE: 21-Feb-2014

2 Summary of Results: HIGH RISK LUMINAL-TYPE (B)

3 MammaPrint Risk-of-Recurrence Result: HIGH RISK

MammaPrint Index: 0.855
High Risk Range: 0.855-1.000

Average 10-Year Risk-of-Recurrence if Untreated: 29%

4 BluePrint Molecular Subtype Result: LUMINAL-TYPE

5 Clinical Pathological High Risk/MammaPrint High Risk*
% of patients without metastasis at 5 years (Distant Metastasis Free Survival)

Endocrine Therapy + Chemotherapy: 90.6%

6 Clinical Pathological Low Risk/MammaPrint High Risk*
% of patients without metastasis at 5 years (Distant Metastasis Free Survival)

Endocrine Therapy Alone: 95.0%
Endocrine Therapy + Chemotherapy: 95.8%

Agendia Inc. Science Park, 1066 Xth Avenue, The Netherlands. T + 31 (0)20 482 1310 customer.service@agendia.com Page 1 of 2

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
Summary of Results

Clinical Risk Assessment in the MINDACT Trial³

ER status	HER2 status	Grade	Node status	Tumor Size	Clinical Risk in MINDACT
ER positive	HER2 negative	Well differentiated (Grade 1)	Node-negative	≤ 1cm	Low
			T+ positive nodes	≤ 1cm	High
		Moderately differentiated (Grade 2)	Node-negative	≤ 1cm	Low
			T+ positive nodes	≤ 1cm	High
		Poorly differentiated or undifferentiated (Grade 3)	Node-negative	≤ 1cm	Low
			T+ positive nodes	≤ 1cm	High
ER positive	HER2 positive	Well differentiated (Grade 1)	Node-negative	≤ 1cm	Low
			T+ positive nodes	≤ 1cm	High
		Moderately differentiated (Grade 2)	Node-negative	Any size	High
			T+ positive nodes	≤ 1cm	Low
		Poorly differentiated or undifferentiated (Grade 3)	Node-negative	≤ 1cm	High
			T+ positive nodes	Any size	High
ER negative	HER2 negative	Well differentiated (Grade 1)	Node-negative	≤ 1cm	Low
			T+ positive nodes	Any size	High
		Moderately differentiated (Grade 2)	Node-negative	≤ 1cm	Low
			T+ positive nodes	Any size	High
		Poorly differentiated or undifferentiated (Grade 3)	Node-negative	≤ 1cm	High
			T+ positive nodes	Any size	High
ER negative	HER2 positive	Well differentiated (Grade 1)	Node-negative	≤ 1cm	Low
			T+ positive nodes	Any size	High
		Moderately differentiated (Grade 2)	Node-negative	≤ 1cm	Low
			T+ positive nodes	Any size	High
		Poorly differentiated or undifferentiated (Grade 3)	Node-negative	≤ 1cm	High
			T+ positive nodes	Any size	High

Survival Stratified by Blueprint Molecular Tumor Subtyping⁴

Breast cancer is a heterogeneous disease and the grouping of breast cancers into distinct clinically-relevant subtypes enables us to inform our treatment decisions on ailing. Blueprint is a functional molecular subtyping assay that classifies breast cancer into three distinct subtypes: Luminal-type, HER2-type and Basal-type by determining the mRNA levels of 80 genes that best discriminate among the following molecular subtypes.^{4,5}

Combining Mammaprint and Blueprint allows patients to be stratified into these subgroups:

- Luminal-type Mammaprint, Low Risk (Luminal A)
- Luminal-type Mammaprint, High Risk (Luminal B)
- HER2-type
- Basal-type

Subtype	Chemotherapy-Free Survival
Low Risk Luminal-type	Low distant rate of pathology complete response (pCR) (6%)
High Risk Luminal-type	High distant rate of pathology complete response (pCR) (10%)
HER2-type	pCR 47%
Basal-type	pCR 37%

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Summary of Results
LATE RECURRENCE (20YR) LOW RISK

MammaPrint Risk-of-Recurrence Result:
LATE RECURRENCE (20YR) LOW RISK

MammaPrint Index: -0.455

MammaPrint Index (MPI)
Patients (Candidates for Endocrine Therapy) (MammaPrint Index < +0.355) | MammaPrint LRLR: Limited Endocrine Therapy Strategy* (MammaPrint Index > +0.355)

100% Predicted Breast Cancer Specific Survival*
Years: 0, 5, 10, 15, 20

11 MPI does not reliably predict for individual therapy. (Late Recurrence (20yr) Low Risk)

Patients with a MPI < +0.355, within the Late Recurrence (20yr) Low Risk range (LRLR), have an exceedingly low likelihood for late recurrence of breast cancer at 20 years* with limited or no endocrine treatment.

Patients with a MPI > +0.355 have shown to benefit from tamoxifen, and are therefore candidates for hormonal therapy.

13 MammaPrint Late Recurrence (20yr) Low Risk^{4,5}

% of patients with Breast Cancer Specific Survival at 20 years⁴

Untreated	94%
Tamoxifen (2-5 years) †	97%

Breast Cancer Specific Survival (BCSS) is the percentage of patients in the study who have not died from breast cancer within the defined period of time (20 years). Patients who died from causes other than breast cancer are not counted in this measurement.

† 65% of patients in the ST03 trial received tamoxifen for any 2 years⁴

Agenda Summary Page Disclaimers:
The Summary Pages are provided for general informational purposes and are not part of any official diagnostic report. Please refer to the official individual patient record for final results. This information (including without limitation, advice and recommendations) and services are neither medical nor health care advice for any individual problem nor a substitute for advice and services from a qualified health care provider familiar with the patient's medical history. All publication information can be found at www.agendia.com.

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Page 1 of 1

PATIENT NAME: **Example Last, First** DOB: 10-Oct-1961 **AGENDIA**

GENDER: Female	ORDERED BY: Physician, Test	REGISTRATION #: Example
SPECIMEN ID: SIG 2	ACCOUNT: Agendia Test -NL- Amsterdam NL	SPECIMEN TYPE: FFPE, Needle Core
PATIENT MRN: MRN 2		SPECIMEN SOURCE: Left Breast
CUSTOMER REF: CRFP 2		COLLECTED DATE: 13-Apr-2020
		RECEIVED DATE: 14-Apr-2020
		REPORTED DATE: 14-Apr-2020

MammaPrint® FFPE Result **High Risk**

The breast cancer tissue sample submitted was analyzed by MammaPrint FFPE, an iVDMA 70-Gene Profile of Breast Cancer for Metastatic Risk that has been validated to correlate with high or low outcome risk for distant metastases in patients with invasive breast cancer. This risk assessment is based on a retrospective analysis of a prospective observational study that included 842 breast cancer patients treated and not treated with adjuvant therapy. Treatment was selected according to clinical assessments that included MammaPrint test results. The risk for distant metastases in unselected patients who did not receive adjuvant treatment was not studied; therefore, MammaPrint FFPE should be used as a prognostic marker only. As a group, "low risk" patients (see those in the MammaPrint FFPE clinical validation (RADIERS) study) have a 1.7% chance (95% CI 0.3-3.1), and "high risk" patients have an 11.7% chance (95% CI 6.6-16.8) that their cancer will recur within 5 years (not accounting for any covariates other than the patient's MammaPrint FFPE status).

Additional Comments:

Assay Description

The U.S. FDA has provided iVDMA clearance of MammaPrint with FFPE tissue for patients with Stage I and II invasive breast cancer, tumor size ≤ 5 cm, lymph node negative, based upon the development and validation of the MammaPrint assay as reported in Nature, New England Journal of Medicine, JNCI, BMC Genomics, Pers. Medicine, and Ann Oncol. The test is performed using a microarray-based gene expression profile that was independently validated on 5-year outcome data on a patient cohort. A FFPE sample's MammaPrint index (MPI) falls within a pre-defined area around the classification cut-off between -0.050 and +0.050, the classification accuracy is less than 50%. See MammaPrint Physician's Brochure found on www.agendia.com for more information.

Sign Off

Jan-Peter Janssen, MD, PhD
 Director of
 In Vitro Diagnostics Use [X]

Agendia Inc. Federal law restricts this device to sale by or on the order of a physician.
 Agendia, Inc. (05D1089250) is certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) as qualified to perform high-complexity clinical testing. MammaPrint FFPE is an IVD with the progress of patients diagnosed with breast cancer. Results of a test should be used in the treatment medical regimen of the testing patient only, and consultation of a medical professional concerning the patient's condition, including other pathologic tests, is essential with the receipt of any test results. MammaPrint test should only be performed using approved samples, using the right, exact equipment, patients to capture the 100% (200 genes) from a gene expression profile. The procedure for sample use is an essential element of the patient sample group. This test is not performed at Agendia's sites. Contact information about MammaPrint FFPE can be found at www.agendia.com.

References:
 1) Drukker CA et al. Nature 2002;415(31):530-536.
 2) Drukker CA et al. Int J Cancer 2013;133(4):929-36.
 3) Van de Vijver MJ et al. New Engl J Med 2002; 347(25):1999-2009.
 4) Buyse M et al. J Natl Cancer Inst 2006; 98(17):1183-1192.
 5) Glas AM et al. BMC Genomics 2006;7:278.
 6) Delahaye LJM et al. Pers Med 2013;10:801.
 7) Mook S et al. Ann Oncol 2010;21(4):717-722.
 8) Mook S et al. Ann Oncol 2010;21(4):717-722.

Agendia NV | Radarweg 60 | 1043 NT Amsterdam | The Netherlands
 Phone: +31 (0) 20 462 1510 | customerservice@agendia.com | www.agendia.com | © 2020

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